



Cap

Mob Cap

Features

Made of soft non-woven material, which gives the wearer a comfortable experience and secures hair and minimizes the risk of contamination. The elastic trim accommodates a variety of hairstyles. There are single elastic string or double elastic strings available.



SKU	Size	Color	Material	Package	Carton Dimension
MC101002	21" single elastic	White	SPP 10GSM	100pcs/pk,10pks/ctn	310*240*250mm
MC101001	21" single elastic	Blue	SPP 10GSM	100pcs/pk,10pks/ctn	310*240*250mm
MC141002	21" double elastic	White	SPP 14GSM	100pcs/pk,10pks/ctn	380*250*260mm
MC141001	21" double elastic	Blue	SPP 14GSM	100pcs/pk,10pks/ctn	380*250*260mm



Bouffant Cap

Features


Made of Spunbonded polypropylene, and sewn with elastic band, with its class one flammability rating and latex free features, it gives your head great protection and comfort. There are size 18",21" and 24" available, which fit most people. It's widely used for cleaning, food service, healthcare, laboratory and homecare.

- Latex free
- Suitable for isolation and basic protection from bacteria and particulate
- Specialized type or design for particular occasions
- Practical barrier for wide range of applications
- Soft and light weight
- Good fit, feel and performance



Applications

 Medical Purpose / Examination

 Industrial purpose / PPE

 Laboratory

 Healthcare and nursing

 General housekeeping

 IT Industry



SKU	Size	Color	Material	Package	Carton Dimension
BC101002	21"	White	SPP 10GSM	100pcs/pk,10pks/ctn	340*270*270mm
BC101001	21"	Blue	SPP 10GSM	100pcs/pk,10pks/ctn	340*270*270mm
BC141002	21"	White	SPP 14GSM	100pcs/pk,10pks/ctn	340*270*270mm
BC141001	21"	Blue	SPP 14GSM	100pcs/pk,10pks/ctn	340*270*270mm

