



Shoe Cover

Features

- Suitable for isolation and basic protection from particulate
- Specialized design for particular occasions
- Practical barrier for wide range of applications
- Soft and light weight
- Good fit,feel and performance



SKU	Size	Color	Material	Package	Carton Dimension
SC351001	160*410mm	Blue	SPP 35GSM	100pcs/pk,10pks/ctn	520*300*320mm
SC401001	160*410mm	Blue	SPP 40GSM	100pcs/pk,10pks/ctn	550*300*320mm

